

## Chlorogenic Acid Analysis

Chlorogenic acid is one of polyphenolic compounds found in coffee.

<b>Samples</b>	Coffee, egg plant, burdock, sweet potato, etc.
<b>Required Sample Amount</b>	30 g
<b>Turnaround Time</b>	2 weeks
<b>Analytical Method</b>	HPLC
<b>Quantification Limit</b>	1 mg/100 g
<b>Standard for Testing</b>	5-Caffeoylquinic acid (CAS No. 327-97-9)

For more information, please contact us at [kaigaisoudan@jfri.or.jp](mailto:kaigaisoudan@jfri.or.jp)