

Curcuminoids Analysis

Curcuminoids are polyphenolic yellow pigments responsible for the yellow color of the root of turmeric.

We can analyze three curcuminoids respectively.

Samples	Turmeric, related processed foods, etc.
Three Curcuminoids	Curcumin, Demethoxycurcumin, Bisdemethoxycurcumin
Required Sample Amount	Liquid: 20 g, Powder: 10 g, Turmeric itself: 100 g
Turnaround Time	3 weeks
Analytical Method	HPLC
Quantification Limit	0.0001 g/100g

For more information, please contact us at kaigaisoudan@jfri.or.jp